



# Teats

Materni's ready to use ETO free range of teats for neonatology, maternity and pediatrics



## Materni safe & easy teat standard product range

Article	Article code	Screw thread	Microbiology	Packaging	Items/Box
<b>SMALL TEAT</b>					
Materni teat - Neonatology / Small - Slow Flow	MEPA0105	EU	CBM	Individually packed	140
<b>MEDIUM TEATS</b>					
Materni teat - Medium - Slow Flow	MEPA0103	EU	CBM	Individually packed	140
Materni teat - Medium - Fast Flow	MEPA0104	EU	CBM	Individually packed	140
Materni teat - Medium - Extra Fast Flow	MEPA0118	EU	CBM	Individually packed	140
<b>LARGE TEATS</b>					
Materni teat - Large - Slow Flow	MEPA0101	EU	CBM	Individually packed	140
Materni teat - Large - Fast Flow	MEPA0102	EU	CBM	Individually packed	140
Materni teat - Large - Extra Fast Flow	MEPA0115	EU	CBM	Individually packed	140

## Materni soft teats **[NEW]**

Article	Article code	Screw thread	Microbiology	Packaging	Items/Box
<b>SOFT TEAT</b>					
Materni Teat - Soft - Large Slow Flow	MXPA0142	EU	CBM	Individually packed	140



Designed fit for function

# Materni teats



Materni's teats are produced and packed, ready to use, in Clean Room ISO Class 5 aseptic sterile conditions. Our Clean-by-Manufacturing (CBM) production and products are now accepted as the new standard by the Belgium and French authorities. We eliminate the use of Ethylene Oxide and gamma irradiation. Our vision: Zero chemicals, local production, sustainable, innovative and affordable.

Materni's teats grow with your baby. There is a suitable size for every baby, for both standard and thicker milk so that your baby can drink easily at all times. The Materni single use teats are a suitable medical device for Oral Rehydration Therapy (ORT).


Finding the perfect teat for each baby is depending on the baby's age, latching skills and the consistency of the preferred milk. Materni offers a broad range of affordable teats with different sizes and speeds, suitable for each situation.

## A total and versatile range of teats for the hospital environment


- CE Class I - Medical device, single use
- European Medical Devices Regulation 2017/745
- Safe and hygienic, individually packed and ready to use
- Zero chemicals philosophy: ETO, BPA, phthalate and latex free.
- Food and medical approved polymers - compliant with regulation (EU) No 10/2011
- Bite resistant compliant according EN 14350 (except small teat)
- No teat collapse possible during feeding because of reinforced design (except small teat)
- Ergonomic: wide teat base
- Fit for different ages: small (premature), medium (< 2.75 kg) large (> 2.75 kg)
- Choose a desired flow function:
  - **Slow Flow (SF)** – hard working – breast feed supportive – thin milk – ORS (Oral Rehydration Solution)
  - **Fast Flow (FF)** – for thicker milk – ORS (Oral Rehydration Solution)
  - **Extra Fast Flow (XFF)** for high viscous formula milk – ORS (Oral Rehydration Solution)
- Fit on bottles with EU DIN screw thread
- Unique one-piece ring & teat design = no leak

## Materni teats have unique advantages

### Opening



### Position



## Small teat

Neonatology (premature)



Slow Flow  
MEPA 0105

## Medium teats

Maternity and pediatrics (< 2.75 kg)



Slow Flow  
MEPA 0103



Fast Flow  
MEPA 0104



Extra Fast Flow  
MEPA 0118

## Large teats

Maternity and pediatrics (> 2.75 kg)



Slow Flow  
MEPA 0101



Fast Flow  
MEPA 0102



Extra Fast Flow  
MEPA 0115

“With such a wide range of teat sizes and flow rates, the Mum and I can choose the teat that best suits the baby's drinking technique.

The slow flow and smaller teats are most useful initially.

Calm feeding is crucial for the digestion and well-being of the baby.”

Experienced midwife Hilde